



# WHEELCHAIR PLATFORM SCALE

## Model: AMW-300



### General Specifications

- Maximum capacity of 300kg.
- 100g resolution provides precise weighing.
- Low profile platform with easy wheelchair access from ramps on both sides.
- Large steel weighing plate and handrail design, makes it ideal to weigh obese patients.
- Easy to read LCD display with white colour LED backlight.
- Lock and Hold functions are available as standard
- Handrails and ramps can easily be removed if necessary.
- Platform size is 880mm x 840mm.
- BMI function to determine underweight/ normal weight / surplus weight.
- Powered by AC adaptor with a rechargeable battery.
- Operating time up to 40 hours without backlight, and the charging time is approximately 12 hours.

#### AMW-300 SCALE SPECIFICATIONS

Capacity	300kg
Resolution	100g

\* Specifications are subject to change without notice