




BARIATRIC STAND-ON SCALE

Model: AMH-400

General Specifications

- Maximum capacity of 250kg or 400kg.
- 100g / 200g resolution provides precise weighing.
- BMI function to determine underweight/ normal weight / surplus weight.
- Easy to read LCD display with white colour LED backlight.
- LCD character height: 25mm.
- Lock and Hold functions are available as standard.
- Platform size is 600mm x 600mm.
- Powered by AC adaptor with a rechargeable battery.
- Operating time up to 40 hours without backlight, and the charging time is approximately 12 hours.



AMH-400 SCALE SPECIFICATIONS		
Capacity	400kg	250kg
Resolution	200g	100g

* Specifications are subject to change without notice