


LOW PROFILE ANIMAL WEIGHING SCALE

Model: AW-300

A & D Weighing offers a robust animal weighing scale perfect for weighing a variety of small animals up to 300kg, including dogs, cats, and other pets. This highly accurate digital scale arrives pre-calibrated, and the mild steel finish allows for easy cleaning and minimal maintenance. This scale base is connected to the AD-4406A indicator display.

General Specifications

- Low profile design with rigid construction
- Dual range capacity 150kg x 50g / 150 - 300kg x 100g
- Portable lightweight body
- Base dimensions are 1100mm x 510mm
- Hold function for animal weighing
- Large liquid crystal display (LCD) with 25mm character height
- Included is a non-slip rubber mat for improved animal stability
- Indicator can be wall or desk mounted with bracket attachment
- Powered by AC adaptor included (or C size batteries, not included)
- Optional RS232C interface



AW-300 SCALE SPECIFICATIONS		
Dual Capacity	0-150kg	150-300kg
Resolution	50g	100g